



The Power of Bonsai® Oxygen Conservation

How does the Bonsai® oxygen conservation savings ratio of up to 6:1 translate into dollar and cents savings for you? The duration chart on the right illustrates how much longer your patient's oxygen supply will last using the revolutionary Bonsai® conserver!

BONSAI®	SETTING	1	2	3	4	5	6	7	CF 2 lpm
Cylinder Type	Cylinder Volume	Estimated Cylinder Duration in Hours (Based on 20 breaths per minute)							
M2	36 liters	2.3	1.4	0.9	0.8	0.7	0.6	0.5	0.3
M4(A)	113 liters	7.2	4.3	2.9	2.4	2.0	1.8	1.6	0.9
M6(B)	164 liters	10.5	6.2	4.3	3.5	3.0	2.6	2.4	1.4
ML6	171 liters	11.0	6.5	4.5	3.7	3.1	2.7	2.5	1.4
M7	198 liters	12.7	7.5	5.2	4.2	3.6	3.2	2.8	1.7
M9(C)	246 liters	15.8	9.3	6.4	5.3	4.5	3.9	3.5	2.1
D	425 liters	27.2	16.1	11.1	9.1	7.7	6.8	6.1	3.5
E	680 liters	43.6	25.8	17.7	14.5	12.3	10.9	9.8	5.7

Estimated up to 6:1 savings dependent on setting. Durations based on cylinder manufacturer's suggested service pressures.

For more information, including more detailed product specifications, contact us at:

01.239.687.1280 • www.chadtherapeutics.com • FAX 01.239.687.1285 • Toll Free 800.423.8870