



NEWS RELEASE

Contact: Len Serafino
Vice President of Sales
Chad Therapeutics
2975 Horseshoe Drive South
Suite 600
Naples, FL 34104
(239) 687-1280, ext. 356

For Immediate Release

CHAD[®] Therapeutics Launches the EVOLUTION[™] Electronic Oxygen Conserver

Naples, Florida -- *CHAD Therapeutics announced today the launch of its new electronic oxygen conserver, the EVOLUTION[™]. The newest in the line of CHAD products, the EVOLUTION[™] is able to provide a minimum of two years operation using just two AA alkaline batteries. With a sensitive triggering mechanism, fast delivery, and conserve settings up to 7 LPM, the EVOLUTION[™] can accommodate a broad range of patients and conditions.*

The single-lumen design provides a 5:1 savings ratio at all settings. Features also include, a uniform volume of oxygen with each pulse ranging from 14-40 breathes per minute and operation between 200 and 3000 PSI.

The functionality and reliability of the EVOLUTION[™] demonstrates CHAD's continued commitment to patient health and satisfaction.

For more information, please contact our customer service department at 1.800.423.8870 or visit www.chadtherapeutics.com.

About CHAD Therapeutics

For over 25 years, CHAD Therapeutics has been the only manufacturer that specializes in oxygen conservation. As the market leader, we offer the most innovative, clinically proven, and complete line of electronic, pneumatic, and disposable oxygen conservers on the market, ensuring you meet the diverse needs and preferences of your oxygen patients.